



I-Can/U-Can

Method	Wheelchair	Lockdown/Suitable
Rear Wheel Drive	Invacare	
	Storm 4	Not Suitable
	Spectra XTR	Q'Straint QLK 100
	Spectra Plus	Q'Straint QLK 100
	Topan	Q'Straint QLK 100
	Stream	Q'Straint QLK 100
	Mirage	Q'Straint QLK 100
Centre Drive	Invacare	
	TDX	Q'Straint QLK 100
	Typhoon 2	Q'Straint QLK 100
	Tornado	Not Suitable
Front or Rear Wheel Drive	Handicare	
	Beatle	New chair set up with reinforced battery box will suit vehicle and will fit with QLK 100
Rear Wheel Drive	Sunrise	
	Salsa	Not Suitable
	Rumba	Q'Straint QLK 100
	F40VR	Q'Straint QLK 100
	Ibis	Not Suitable
Centre Drive	Sunrise	
	Groove	Q'Straint QLK 100
Centre Drive	Pride	
	Q600 and Q610	Dock and Lock
	Jazzy 1121 and 1103	These need to be checked
Centre Drive	ETAC	
	Balder	Balder Spigot
Rear Wheel Drive	Auto Adapt	
	Carony Drive	Carony Drive Lock
	Permobil	
	C300 and C400 and C500	These need to be checked
Centre Drive	Pride	
	600 and 610	Q'Straint QLK 100

ALWAYS REMEMBER THAT THE HEIGHT OF THE CUSTOMER MUST BE TAKEN INTO CONSIDERATION